



Part #: BSN-150



1-Year Sensor Pad  
(included)



30-Day Sensor Pad  
(optional)

## Economy Bed Monitoring System

Curbell proudly introduces our new **Economy Bed Monitoring System**, a low cost alternative to our electrically powered Bed Monitoring System. It has similar styling to our Chair Monitor, with the added benefit of a flashing LED indicator that indicates the unit is on or off. If the resident attempts an unsafe exit from the bed, the monitor sounds an audible alarm. The unit is portable from bed to bed and is a fall management product that offers exceptional value.

The **Economy Bed Monitoring System** comes with a 1-Year Reusable Pad or if desired, a Disposable 30-Day pad is also available (sold separately).



7 Cobham Drive • Orchard Park, NY 14127-4180  
Toll-Free: (800) 235-7500 • Fax: (716) 667-7775  
Web: [www.fallmanagement.com](http://www.fallmanagement.com)  
Email: [electronics@curbell.com](mailto:electronics@curbell.com)



## STANDARD FEATURES

- Battery operated
- Portable from bed to bed
- Visual LED indicator
- 1-Year Reusable bed sensor pad
- Adjustable volume control
- Several mounting options (Velcro®, clip & wire holder)

*Note: Disposable 30-Day pads are also available*

## TECHNICAL SPECIFICATIONS

- Power supply:  
Standard 9-volt lithium battery
- Control unit dimensions:  
4.38" deep x 3.25" wide x 1.5 wide"
- Weight:  
approximately 6 oz.
- Mounting:  
Velcro®, clip or wire holder.
- Audio Alarm: Intermittent audible alarm tone @ 98db, central frequency of 2400+/-500 Hz.
- Low battery indicator:  
Intermittent chirp every 10 seconds
- Sensor dimensions:  
3.5" wide x 29" long

## LIMITED WARRANTY

- 2-Year on Control Box
- 1-Year on Reusable Pads
- 30-Days on Disposable Pads