

Chair Monitoring System
Part #: BSN-065

Product Specifications

The **Chair Monitoring System** is a battery-operated, portable monitor with an audible exit alarm. To use, the reusable 1-year chair sensor pad is placed on the seat of the chair. The resident is then seated and the monitor turned on for use. If the resident attempts to leave the chair, the monitor sounds an audible alarm.

Complete System Components:

Control Box, 1-year reusable bed sensor pad and lithium battery



Standard Features:

- Battery Operated
- Three mounting options: (Velcro, clip and wire hanger)
- Low battery warning
- In-room alarm on / off
- Adjustable alarm tone volume
- Automatic reset

Other Important Features:

- Made in USA
- Portable
- 2-year control box limited warranty
- 1-year reusable sensor pad limited warranty

Technical Specifications:

Control Unit Dimensions	4.38" deep x 3.25" wide x 1.5" long
Weight	Approximately 6 oz
Power Source	Standard 9-volt lithium battery
Audio	Audible Alarm tone
Chair Sensor Pads	Reusable (1-year)